

Labmex™
INNOVATIVE DESIGN FOR LABORATORY



SafeMex

Laboratory Fume Cupboard
for General Purpose

FEATURES

SAFETY

Labmex's SafeMex fume cupboard is tested and comply with European Standard EN14175-3

Sash lock to restrict user to work at maximum safety sash position

Optional flow sensor and display for additional air flow safety monitoring

STRUCTURE

- Epoxy coated steel coating for chemical resistance
- Baffle systems for low turbulence inflow
- Phenolic resin, Epoxy resin or Ceramic worktop available
- Base cabinet with key lock option

EXTERNAL DIMENSION

	Item No.	Dimension (W x D x H mm)
4 FT	LX-SFM-1200	1200 x 895 x 2355
5 FT	LX-SFM-1500	1500 x 895 x 2355
6 FT	LX-SFM-1800	1800 x 895 x 2355