



LABSIM

Static laboratory bench

FEATURES

- 20mm phenolic resin worktop with anti-microbial surface are chemically, scratch and abrasion resistant
- Work Top edges are polished with no sharp corner
- Rectangle Frame with 40mm x 60mm x 2mm welded box section construction
- Oven baked powder coated finish steel
- 30mm Adjustable levelling feet
- Easy to Install with just one tool to fix
- Loading capacity >100kg



SPECIFICATIONS

Model No.	Dimension WxDxH (mm)
LX-LS-1000	1000x750x844
LX-LS-1200	1200x750x844
LX-LS-1500	1500x750x844
LX-LS-1800	1800x750x844