



LABSIM

Laboratory Seating

FEATURES

- Materials meets BS7176:2007 + A1:2011 for resistance to ignition of upholstered furniture
- Chemical resistant seating material for daily cleaning with disinfectant
- Gas lift height adjustment for comfort sitting and ergonomic
- Smooth castors allows to move smoothly without rub against the flooring
- Various choice of designs that fit different purpose



The All-rounder

The Core

The Fundamental

The Pose

The Pony

Our Products

		Ordering code
The Fundamental	Best selling upholstered vinyl backrest chair with adjustable seat height 610 – 800mm, ergonomic mechanism allowing backrest adjustment, including height, depth and tilt angle. Come with anti-slip design plastic foot ring for extra comfort.	LX-LS-VBC-S-PR
The All-rounder	A traditional upholstered vinyl stool, versatile for different situation and purpose. Adjustable seat height 600 – 720mm, with metal foot ring for extra comfort.	LX-LS-VST-T-MR
The Core	A easy to clean and durable polyurethane backrest chair with adjustable seat height 610 – 800mm, ergonomic mechanism allowing backrest adjustment, including height, depth and tilt angle. Come with metal foot ring for extra comfort.	LX-LS-PBC-T-MR
The Pose	Polyurethane back support stool design for laboratory usage, it is space saving but do not sacrifice comfort. Adjustable seat height 620 – 730mm, with metal foot ring for extra comfort.	LX-LS-PLB-T-MR
The Pony	Upholstered vinyl saddle design stool allowing sitting with a natural curve, it is specifically designed to ease back problems and strengthens back muscles. Adjustable seat height 510 – 690mm, with metal foot ring for extra comfort.	LX-LS-VSD-S-MR

Available options: lower/higher seat height range, wheels/foot glide and foot ring material