



FabMex

Laboratory Fume Cupboard
in Fiberglass Reinforced Polyester (FRP)



FEATURES

FRP is a strong structural strength, superior in chemical and heat resistance material. It has good resistance to most solvents, bases and acids and high fire resistance.

Our Fume Cupboard constructed molded FRP, and it have baffles for smooth air exhaust with low turbulence.

It has a make up air design enhance good ventilation when there is insufficient air provided.

Sash glass is 5mm tempered glass with PP made frame and sash handle

Worktop Material:

Choose from Phenolic resin, Epoxy Resin and Ceramic worktop

EXTERNAL DIMENSION

	Item No.	Dimension (W x D x H mm)
4 FT	LX-FFM-1200	1245 x 900 x 2300
5 FT	LX-FFM-1500	1500 x 900 x 2300
6 FT	LX-FFM-1800	1800 x 900 x 2300